


TESTICULAR IMPLANTS

SMOOTH




Benefits

- Egg-shaped to closely match the morphology of a testicle, prefilled with Naturgel 1.
- A range featuring five sizes, smooth surface, corresponding to the patient's morphology.
- A double occlusion patch, containing an eye to help suturing the implant.
- A removable thread to guide the needle and prevent a puncture of the implant.



Height (mm)	Diam. (mm)	Vol. (ml)	Reference
25	21	6	LS 01 025
34	27	12	LS 01 034
40	30	17	LS 01 040
50	33	27	LS 01 050
53	38	38	LS 01 053



Our sterile, single-use testicular implants are sold by the unit. They are class IIb devices.

CALF IMPLANTS

SMOOTH



Benefits

- Symmetrically-shaped to create a natural curve on the inner or outer side of the calf.
- Prefilled with Naturgel 1.5 and a smooth surface providing a natural feeling.
- Provided with an insertion bag to ease introduction.



Length (mm)	Width (mm)	Proj. (mm)	Vol. (ml)	Reference
120	40	15	40	LS 03 120
175	35	20	60	LS 03 175
180	42	21	85	LS 03 180
185	46	23	110	LS 03 185
190	50	25	140	LS 03 190
195	57	27	165	LS 03 195



Our sterile, single-use calf implants are sold by the unit. They are class IIb devices.