

# ROUND, BICONVEX ) & ANATOMICAL GLUTEAL IMPLANTS



SEBBIN PARIS  
ESTHETIQUE & RECONSTRUCTION

Recreating harmonious bodies for a new vision of rebirth.

# Anatomical gluteal implants




Reference	Volume (mL)	Width (mm)	Height (mm)	Projection (mm)
LS 05 185	185	97	139	30
LS 05 260	260	107	148	32
LS 05 335	335	117	162	36
LS 05 400	400	129	174	42

This shape is recommended for patients when:

- The majority of the gluteal volume is in the upper part of the buttock,
- The gluteal muscle is long,
- The buttocks form a high square,
- The patients are tall.

It mainly allows to increase the volume of the lower part of the buttock.




# Round gluteal implants

Reference	Volume (mL)	Base (mm)	Projection (mm)
LS 04 200	200	108	30
LS 04 235	235	113	33
LS 04 270	270	120	34
LS 04 300	300	125	36
LS 04 330	330	127	38
LS 04 360	360	128	42
LS 04 400	400	132	43
LS 04 440	440	136	44

This implant is recommended in the following indications:

- If the majority of the gluteal volume is in the lower part of the buttock,
- If the gluteal muscle is short,
- If the buttocks form a small square.




# Biconvex gluteal implants



Reference	Volume (mL)	Base (mm)	Projection (mm)
LS 06 370	370	127	52
LS 06 410	410	130	53
LS 06 480	480	136	56
LS 06 530	530	143	59

- The biconvex shape answers the same indications as the round shape but ensures more projection for an identical base.
- Moreover, in the event of a shift from deep face to superficial face, the biconvex implant being in all points symmetrical, makes this displacement indiscernible for the patient.



# Features



Groupe SEBBIN is one of the first European manufacturers to propose this implantable medical device.

The gluteal implants prefilled with silicone gel are more flexible, easier to introduce and offer a better projection compared to the silicone elastomer monobloc implants.

During the introduction, the implant can be strongly stressed; thanks to their specific very high cohesive silicone gel, our gluteal implants return rapidly to their initial shape.

Our sterile gluteal implants, for single use, are sold by unit. They can be manufactured upon request with a textured envelope under the reference LST.

The dimensions of our round, biconvex and anatomical implants can differ slightly from the data mentioned in this brochure.



# Our commitments



## CONTINUITY

Production know-how mastered for almost 30 years.

## RELIABILITY

Biological and mechanical evaluation of end products.  
Visual controls on 100 % of implants at each step of the manufacturing process.

## TRACKING

Each implant is identified by a unique serial number allowing tracking from the raw materials used to its implantation.

## PEACE OF MIND

Highly qualified personnel working in clean room.

This document is dedicated to health care professionals. According to the Medical Device Directive 93/42/EEC, the gluteal implants are Class IIb medical devices manufactured by Groupe SEBBIN and are intended to be used in the context of a plastic, reconstructive and cosmetic surgery. These devices are CE marked by the notified body MDC number 0483. Please read carefully the instructions for use before using. These devices may not be reimbursed by the health insurance organizations.

### MANUFACTURER:



**GROUPE SEBBIN SAS**  
39, Parc d'Activités des Quatre Vents  
95650 Boissy l'Aillierie - France  
Tel.: + 33 1 34 42 13 28  
Fax: + 33 1 34 42 16 88  
[www.sebbin.com](http://www.sebbin.com)

### INTERNATIONAL OPERATIONS:

**Sebbin Benelux:**  
Tel.: +32 3 451 26 81 - Fax: +32 3 451 26 89  
Email: [info@laboratoires-sebbin.be](mailto:info@laboratoires-sebbin.be)  
**Sebbin Deutschland GmbH:**  
Tel.: +49 2102 70 88 00 0 - Fax: +49 2102 70 88 00 99  
Email: [service@sebbin-implantate.de](mailto:service@sebbin-implantate.de)  
**Sebbin France:**  
Tel.: +33 1 34 42 13 28 - Fax: +33 1 34 42 16 88  
Email: [serviceclients@laboratoires-sebbin.fr](mailto:serviceclients@laboratoires-sebbin.fr)  
**Sebbin Iberica Productos Medicos, S.L.:**  
Tel.: +34 93 159 3003 - Fax: +34 93 676 3053  
Email: [info@sebbin.es](mailto:info@sebbin.es)