

INFLATABLE MAMMARY IMPLANTS



SEBBIN PARIS
ESTHETIQUE & RECONSTRUCTION

Recreating harmonious bodies for a new vision of rebirth.





Dr Jenny's diaphragm valve

Recommended Volume (ml)	Base (mm)	Projection (mm)	Reference	
			Smooth	Textured
70	74	26	LS 20 070	LS 21 070
90	90	28	LS 20 090	LS 21 090
120	97	29	LS 20 120	LS 21 120
140	103	32	LS 20 140	LS 21 140
180	113	34	LS 20 180	LS 21 180
220	120	36	LS 20 220	LS 21 220
240	124	38	LS 20 240	LS 21 240
280	127	40	LS 20 280	LS 21 280
300	130	41	LS 20 300	LS 21 300
360	137	43	LS 20 360	LS 21 360
380	141	47	LS 20 380	LS 21 380
430	147	49	LS 20 430	LS 21 430
520	149	55	LS 20 520	LS 21 520
620	150	61	LS 20 620	LS 21 620
720	152	74	LS 20 720	LS 21 720

Not every device will conform to the dimensions given. Slight variations may occur. Only authorized content to be used : sterile saline solution at 0,9% NaCl.

- Large base for natural result
- Wide choice of volumes
- Aesthetic surgery : augmentation
- Reparative surgery : hypoplasia, asymmetry and reconstruction



Our sterile and single use inflatable mammary implants are sold by unit with their filling tube

The filling tube, with the reference of LS 99 100, can also be supplied separately.

CONTINUITY Manufacturing process developed throughout the last 25 years.

RELIABILITY Chemical and biological evaluation of raw materials then microbiological evaluation of each batch. Optical and technical controls of the implants at each step of the manufacturing process.

TRACEABILITY Each implant identification number allows traceability from the raw materials up to the implantation. We introduce our new bar code system GSI-128 to improve everyday life.

SERENITY Highly qualified staff working in clean room.



PVP-SAS 39, Parc d'Activités des Quatre Vents - 95650 Boissy l'Aillierie - France
 Tél. : + 33 1 34 42 13 28 - Fax : + 33 1 34 42 16 88 - www.sebbin.com
 CE 0483 in accordance with Medical Device Directive 93/42 EEC.